

3 simple steps towards diabetes control

1. Test



Insert a test strip into the bolus advisor in the direction of the arrow. The system will turn on automatically.



The 'Apply Sample' screen appears. The bolus advisor is now ready for monitoring.



Use the Accu-Chek Multiclix lancing device to get a blood sample (see overleaf).



Touch the blood drop to the front edge of the yellow window of the test strip. When the test strip has enough blood, the 'Analyzing' screen appears.



The result appears on the display. Approximately 3 seconds later, the detailed 'bG Result' screen appears. Remove and discard the used test strip.

2. Enter Carbohydrates



On the 'bG Result' screen you can add information about the Meal Time, planned food intake or your state of health.


To enter Carbs:

Select the Carbs entry field and press the 'enter' button (shown).

Set the number of carbs and press the 'enter' button again.




3. Receive bolus advice



Press button shown to get to the 'Bolus Advice' screen.

Confirm the entries by pressing 'Confirm'.



The bolus amount is displayed on the screen. Inject the recommended bolus using an insulin pen or syringe.

Select 'OK' by pressing button shown to go to the Main Menu.



When you open a new vial of strips



Switch off the meter & if present remove and discard white chip. Insert the black activation chip from the vial.



Check the use by date on the vial before use. Do not use test strips past the use by date. Once installed the black chip does not need to be replaced.