



Fingers love Accu-Chek® lancets.




Make sure you change to a new lancet for every blood sample.




Changing the lancet every time you take a blood sample is important, the lancet penetrates the skin and to maintain hygiene it is recommended that a fresh lancet is used each time you take a blood sample. Every new lancet in the drum is sterile and after use can be discarded without being touched. The ② on the lancet pack indicates that lancets are for 'single use only'. The same device must not be shared or used between multiple people.

Accu-Chek® Multiclix uses gentle Clixmotion® technology



Clixmotion technology

- No vibration - neither in nor out
- Lancet is brought to a gentle stop and immediately retracted



Conventional technology

- Vibration causes painful tissue damage
- Lancet stops abruptly and bounces back

Experience what's possible.

www.accu-check.com.au - Accu-Chek Enquiry Line 1800 251 816

For people with diabetes. Use only as directed. Consult your healthcare professional for advice. Accu-Chek lancing devices are for single patient use only. Roche Diagnostics Australia Pty Ltd., 31 Victoria Avenue, Castle Hill NSW 2154 ABN 29 003 001 205. ACCU-CHEK, ACCU-CHEK MULTICLIX and CLIXMOTION are trademarks of Roche. All other trademarks are the property of the respective owners. © 2012 Roche Diagnostics.



Fingers love Accu-Chek® lancets.



The Accu-Chek® Multiclix system.



- 1 Press priming button down like a ballpoint pen. The transparent release button turns yellow.



- 2 Press the lancing device firmly against side of fingertip and press release the button.




- 3 To advance to the next lancet, turn priming button as far as it will go and turn back again. The number on the display indicates how many unused lancets are left in the drum.






Replacing lancet drum.



- A Loading a new drum**
Remove cap from the end of your Multiclix device.



- B** Insert a new drum with the blue ring end first. Push to snap the drum into place. Replace cap.



- C Exchanging a used drum**
When number 1 appears in the display, the last lancet needle is in use. When used, remove the cap from the end of your Multiclix device.



- D** Remove used drum. *(It cannot be reinserted.)* You can dispose of the old drum in household waste. Follow Step 'A' to load a new drum.

Changing depth setting.



- Set penetration depth**
Turn the ring near the cap to select depth setting.

At first, select a shallow depth, gradually increasing it until sufficient blood is obtained.
(Start with level 2.)


You're ready to test: **Press, click, and twist.**

ACCUCHEK®



Tips for virtually pain free¹ lancing.

Lancing correctly helps to keep discomfort to a minimum. You can help your patients monitor more comfortably by explaining these five easy steps:



Ensure hands are **washed and clean.**




Lance on the **side of the fingertip** rather than the pad.



Keep the skin taut by pressing the lancing device firmly against the fingertip.



Select a penetration depth set as shallow as possible but on a depth that still produces blood.



Alternate fingers **daily** and take steps to ensure **good blood circulation.**

¹Jendrike N et al. Pain sensation at Fingertips and Palm by using different Blood Glucose Monitoring Systems. Diabetes Technology Meeting. San Francisco, 5-7 November 2009.

ACCUCHEK®